On-X[®]

Instructions for Use CHORD-X ePTFE SUTURE

Nonabsorbable Monofilament
for Chordae Tendineae repair or replacement

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CE
0459

 *On-X* life technologies, inc.[™]
designed for life

TABLE OF CONTENTS

1. DESCRIPTION.....	3
2. INDICATIONS.....	3
3. CONTRAINDICATIONS	3
4. WARNINGS	3
5. PRECAUTIONS	4
6. RECOMMENDED TECHNIQUES	6
6.1 Suture Placement:.....	6
6.2 Establishing Chordae Tendineae Length:.....	7
6.2.1. References:.....	8
7. ADVERSE REACTIONS	8
8. HOW SUPPLIED.....	8
9. DEFINITIONS.....	9
FIGURES	78

INSTRUCTIONS FOR USE CHORD-X EPTFE SUTURE

Nonabsorbable Monofilament for Chordae Tendineae repair or replacement

1. DESCRIPTION

Chord-X ePTFE Suture for Chordae Tendineae Repair or Replacement is a nonabsorbable, monofilament suture manufactured from expanded polytetrafluoroethylene (ePTFE) that has been expanded to produce a porous microstructure which is approximately 50% air by volume. The suture is undyed and contains no additives.

2. INDICATIONS

Chord-X ePTFE Suture for Chordae Tendineae Repair or Replacement is indicated for the repair or replacement of chordae tendineae.

3. CONTRAINDICATIONS

This device is contraindicated for use in ophthalmic surgery, microsurgery, and peripheral neural tissue.

4. WARNINGS

The replacement of chordae tendineae requires accurate sizing to ensure satisfactory functioning of the valve mechanism post-

repair or replacement. Chordae tendineae replaced with sutures that are either too long or too short can result in failure to correct valvular incompetence. See the Recommended Techniques section for determining correct suture length.

Do not use in the presence of infectious endocarditis, which can result in failure of the valve repair.

Device has been sterilized using ethylene oxide. Do not re-sterilize.

Potential tissue invasion of the Chord-X ePTFE Suture may result in the attachment of the suture to the tissue it penetrates. Such attachment could make removal of the Chord-X ePTFE Suture difficult.

5. PRECAUTIONS

Successful use of this device requires appropriate training and experience in cardiac surgical techniques including mitral valve repair.

Do not alter the usual practice of pre-, peri-, or post-operative administration of anticoagulation therapy, which could increase the risk of adverse events associated with this therapy.

Misuse of this suture, like any other suture, can result in severe injury or death to the patient. As with

any suture, care should be taken to avoid damage when handling. Avoid crushing or crimping the suture with surgical instruments or exposing the suture to sharp edges, which may adversely affect the tensile strength of the suture.

At least seven (7), equally tensioned, flat square throws are required to produce a secure knot when tying Chord-X ePTFE Suture for Chordae Tendineae Repair or Replacement. Additional throws may be necessary depending on surgical circumstances. When tying knots with Chord-X ePTFE Suture, tension should be applied by pulling each strand of the suture in opposite directions with equal force. As the knot is tensioned, the air in the suture is forced out. Care should be taken to avoid using a jerking motion which could break the suture. Uneven tensioning of a well-formed square knot may result in an unsecure knot. When the Chord-X ePTFE Suture is properly tensioned and formed, standard surgical knotting techniques will produce a secure knot.

The Chord-X ePTFE Suture device is designed for single use only; do not reuse device. On-X Life Technologies does not have data regarding reuse of this device.

Reuse may cause device failure or procedural complications including device damage, compromised device biocompatibility, and device contamination. Reuse may result in infection, serious injury, or patient death.

6. RECOMMENDED TECHNIQUES

6.1 Suture Placement:

1. After determining the number of chordae tendineae that will be replaced, a Chord-X ePTFE Suture is placed through the papillary muscle, with a pledget placed on each side of the muscle (Figure 1). The suture is then tied. This knot should be tied loosely to avoid ischemia to the tip of the papillary muscle.
2. After all of the sutures have been anchored to the papillary muscle, the two sides of each suture are run through the free margin of the leaflet (approximately 2mm apart). It is useful to increase the friction of the suture in the leaflet to minimize suture slippage during the suture length adjustment steps. Therefore, pass each side of the suture through the leaflet twice (Figure 2). An autogenous or thin ePTFE pledget (not supplied) may be placed on the atrial side of the leaflet for reinforcement (Figure 3).

6.2 Establishing Chordae

Tendineae Length:

1. Once the remaining valve repair procedures are completed, and the length of each chordae tendineae replacement suture is then established. Suture lengths are set to maintain the opposing points of the two leaflets at the same level in the ventricle.
2. Saline solution is injected into the ventricular cavity to close the valve (with aortic root venting for de-airing). If residual regurgitant jets are present or the closure line is not symmetrical, adjust the length of the Chord-X Sutures.
3. Knot the Chord-X ePTFE Suture ends using at least seven (7) flat square throws. Care should be taken to avoid altering the length of the suture while knotting. Shortening of the artificial chordae while knotting could result in over-correction of the leaflet prolapse and restrict the leaflet motion.
4. After the patient is off cardiopulmonary bypass, transesophageal echocardiography should be used to confirm the results of the hydrodynamic testing and evaluate the valve repair.

6.2.1. References:

Zussa, Claudio, Different applications of ePTFE valve chordae: surgical techniques. *J Heart Valve Disease*, 1996; 5:356-361.

Frater, RWM, Gabbay, S, Shore, D, et al. Reproducible replacement of elongated or ruptured mitral valve chordae. *Ann Thorac Surg*, 1983; 35; 14-28.

Zussa, Claudio, *Artificial Chordae in Mitral Valve Surgery*. Austin: R.G. Landes Company, 1994: 159.

7. ADVERSE REACTIONS

Potential adverse effects associated with the use of any suture include: wound dehiscence, infection, and localized transitory inflammatory tissue reaction.

8. HOW SUPPLIED

Chord-X ePTFE Sutures for Chordae Tendineae Repair or Replacement are available as sterile strands in USP 2-0 and USP 3-0 with either 1/2 or 3/8 circular taper point needles and without pledgets.

Chord-X ePTFE Sutures for Chordae Tendineae Repair or Replacement are supplied **STERILE** unless the integrity of the package has been compromised. Do not use the

device if the package has been compromised. The device is for single use only.

9. DEFINITIONS

Use by

Attention, See instructions for use

Do Not Reuse

Catalogue Number

Batch Code

European Authorized Representative

Contents sterile unless package has been opened

Contents sterile unless enclosed package has been opened or damaged. Sterilized by ethylene oxide

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On-X Life Technologies, Inc.
Austin, Texas 78752 U.S.A.