
On-X® Chord-X™ Pre-Measured Loops

For Mitral Chordal Replacement

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On-X® Chord-X™ Suture System

For Mitral Chordal Replacement

ePTFE Nonabsorbable Monofilament Suture, U.S.P. with PTFE Pledgets

1. DESCRIPTION

Chord-X™ Pre-Measured Loops and Chord-X™ Suture System for Mitral Chordal Replacement is an artificial chordal replacement prosthesis configured from nonabsorbable, monofilament suture manufactured from expanded polytetrafluoroethylene (ePTFE). The ePTFE Suture has been expanded to produce a porous microstructure which is approximately 50% air by volume. The suture is undyed and contains no additives. The Chord-X™ Pre-Measured Loops and the Chord-X™ Suture System for Mitral Chordal Replacement prostheses have various product configurations comprised of different suture diameter, loop length and needle configurations.

(See How Supplied)

The Chord-X™ Pre-Measured Loops and Chord-X™ Suture System prostheses are offered: (1) configured with pre-measured, pre-knotted ePTFE suture loops **(See Figure 1)** and (2) configured with ePTFE suture strands with no pre-tied knots allowing the surgeon to adjust and create artificial loops of a customizable length. **(See Figure 2)**

Both the Chord-X™ Pre-Measured Loops and the Chord-X™ Suture System prostheses come equipped with double-armed taper point needles attached. Each prosthesis has 2 PTFE pledgets for papillary attachment. The Chord-X™ Pre-Measured Loops prosthesis features 3 chordal suture loops for the leaflet, complete with needles and pre-tied knots, and 1 suture pair with an attached pledget for the papillary muscle. The Chord-X™ Suture System prosthesis version features 3 pairs of chordal suture strands for the leaflet and 1 pair of suture strands with an attached pledget for the papillary muscle. Each prosthesis is positioned on a series of base-plates for suture management. The base-plates protect the suture from getting kinked, bent or otherwise damaged. If fewer chordal sutures are desired, excess segments can be removed without jeopardizing the rest of the prosthesis.

2. INDICATIONS

Chord-X™ Pre-Measured Loops and Chord-X™ Suture System for Mitral Chordal Replacement are indicated for the repair or replacement of chordae tendinae.

3. CONTRAINDICATIONS

This device is contraindicated for use in ophthalmic surgery, microsurgery, and peripheral neural tissue.

4. WARNINGS

The replacement of chordae tendinae requires accurate sizing to ensure satisfactory functioning of the valve mechanism post repair or replacement. Chordae tendinae replaced with sutures that are either too long or too short can result in failure to correct valvular incompetence. See the **Recommended Techniques** section for determining correct length.

Caution and surgical judgment is necessary when considering the use of this product in the setting of infectious endocarditis. The surgeon must be aware that prosthetic material may become infected in the presence of infectious endocarditis. Use in the presence of infectious endocarditis can result in failure of the valve repair.

The prosthesis has been sterilized using ethylene oxide. Do not re-sterilize.

Potential tissue in-growth of the Chord-X™ Pre-Measured Loops or the Chord-X™ Suture System for Mitral Chordal Replacement made of ePTFE Nonabsorbable Monofilament suture and/or the pledgets made of PTFE felt may result in the attachment of the suture to the tissue it penetrates. Such attachment could make removal of the Chord-X™ Pre-Measured Loops or Chord-X™ Suture System sutures made of ePTFE Nonabsorbable Monofilament suture difficult.

5. PRECAUTIONS

Successful use of this device requires appropriate training and experience in cardiac surgical techniques including mitral valve repair.

Do not alter the usual practice of pre-, peri-, or post-operative administration of anticoagulation therapy, which could increase the risk of adverse events associated with this therapy.

Misuse of this suture prosthesis, like any other suture, can result in severe injury or death to the patient. As with any suture, care should be taken to avoid damage when handling. Avoid crushing or crimping the suture with surgical instruments or exposing the suture to sharp edges, which may adversely affect the tensile strength of the suture. Avoid using knot pushers with imperfections such as burrs or flashing which may damage the suture. Avoid untying the pre-tied suture knots on each pre-measured loop.

At least seven (7), equally tensioned, flat square throws are required to produce a secure knot when tying the Chord-X™ Pre-Measured Loops or the Chord-X™ Suture System for Mitral Chordal Replacement made of ePTFE suture. Additional throws may be necessary depending on surgical circumstances. When tying knots with the Chord-X™ Pre-Measured Loops or the Chord-X™ Suture System prostheses made of ePTFE suture, tension should be applied by pulling each strand of the suture in opposite directions with equal force. As the knot is tensioned, the air in the suture is forced out. Care should be taken to avoid using a jerking motion which could break the suture. Uneven tensioning of a well-formed square knot may result in an unsecure knot. When the Chord-X™ Pre-Measured Loops or the Chord-X™ Suture System for Mitral Chordal Replacement is properly tensioned and formed, standard surgical knotting techniques will produce a secure knot.

The Chord-X™ Pre-Measured Loops and the Chord-X™ Suture System prostheses are designed for single use only; do not reuse device. On-X Life Technologies, Inc. does not have data regarding reuse of this device.

Reuse may cause device failure or procedural complications including device damage, compromised device biocompatibility, and device contamination.

Reuse may result in infection, serious injury, or patient death.

6. RECOMMENDED TECHNIQUES

6.1. Determine the artificial Chord-X™ Pre-Measured Loops or Chord-X™ Suture System Adjustable loop length.

Determine the artificial Chord-X™ Pre-Measured Loops or Chord-X Suture System adjustable loop length via surgeon preferred technique, to maintain the coapting surfaces of the two leaflets at the same level in the ventricle.

6.2. Unpack Chord-X™ Pre-Measured Loops or the Chord-X™ Suture System prostheses from the base-plate.

After determining the artificial chordae loop length, choose the appropriate Chord-X™ Pre-Measured Loops prosthesis or the Chord-X™ Suture System. Open the package using sterile technique. Each prosthesis is positioned on a series of base-plates for suture management. If fewer chordal sutures are desired, excess segments can be removed without jeopardizing the rest of the prosthesis.

1. Remove foam base-plates from case.
2. Unwrap the papillary suture with attached pledget from the baseplate.
3. Unwrap chordal prosthesis sutures 1, 2 and 3 from their respective base-plates.
4. Remove the additional PTFE pledget from case if desired.

6.3. Removing extra loops.

The Chord-X™ Pre-Measured Loops prosthesis version and the Chord-X™ Suture System prosthesis version feature 1 suture pair for a papillary muscle and 3 independent pre-measured chordal suture loops or 3 suture strands (for the leaflet) complete with needles. If fewer chordal sutures or strands are desired, excess segments are easily removed without jeopardizing the rest of the prosthesis.

1. Removing extra loops from the Chord-X™ Pre-Measured Loops prosthesis version:

Remove undesired prosthetic loops by cutting a **single** strand of suture between the knot and the papillary suture strand and discard.

(See Figure 3)

Caution: Be aware if 2 strands are cut, there will be a free-floating piece of suture that will need to be retrieved.

2. Removing extra suture strands from the Chord-X™ Suture System Adjustable prosthesis version:

Remove undesired suture strand by cutting a single chordal suture strand and discard. (See Figure 4)

Caution: Be aware if 2 strands are cut, there will be a free-floating piece of suture that will need to be retrieved.

6.4. Papillary Suture Placement:

Place the Chord-X™ Pre-Measured Loops or the Chord-X™ Suture System prosthesis papillary suture with attached needle through the papillary muscle, with a pledget placed on each side of the muscle. One pledget is pre-attached to the prosthesis. The second pledget is provided in the case. The papillary suture is then tied using at least seven (7) flat square throws. Additional throws may be necessary depending on surgical circumstances. (See Figure 5)

6.5. Leaflet Loop Placement:

1. Chord-X™ Pre-Measured Loops prosthesis version leaflet placement and knot tying:

Once the suture has been anchored to the papillary muscle, pass each required chordal pre-measured suture loop through the free edge of the prolapsed leaflet, taking care to not entrap native chordae as the Chord-X™ Pre-Measured Loops are individually positioned and tied.

The two suture strands of each Chord-X™ chordal Loops are run through the free margin of the leaflet, then tied.

Tie the Chord-X™ Pre-Measured Loop ends using at least seven (7) flat square throws. Continue placement of the remaining Chord-X™ Pre-Measured Loops as needed to complete artificial chordal replacement, tying each one as described. (See Figure 6)

2. Chord-X™ Suture System Adjustable prosthesis version leaflet placement and knot tying:

Once the suture has been anchored to the papillary muscle, pass the required chordal suture strands through the free edge of the prolapsing leaflet, taking care to not entrap normal chordae as the Chord-X™ loops are individually positioned and tied. Establish the length of each chordae tendinae replacement suture via surgeon preferred method. Suture is set to maintain the coapting surfaces of the two leaflets at the same level in the ventricle. **(See Figure 7)**

Surgeons may employ a variety of intraoperative tests to assess chordal length before closing the left atrium. In one such test, the ventricle is filled with fluid under pressure by insufflating it with saline, giving antegrade cardioplegia, or performing both maneuvers. With such testing, if residual regurgitant jets are identified, the surgeon may adjust the length of the Chord-X™ Suture System suture strands.

Knot the Chord-X™ Suture System prosthesis suture ends using at least seven (7) flat square throws. Continue placement of the remaining Chord-X™ Suture System sutures strands as needed to complete artificial chordal replacement, tying each one as described.

Shortening of the artificial chordae while knotting could result in overcorrection of the prolapsed leaflet and restrict the leaflet motion.

6.6. Final Procedures:

Utilize surgeon preferred method or imaging to confirm proper loop length and proper coaptation of leaflets and successful valve repair. After the patient is off cardiopulmonary bypass, transesophageal echocardiography should be used to confirm the results of the hydrodynamic testing and evaluate the valve repair.

Finished Chordal Replacement. **(See Figure 8)**

6.7. References:

- David, TE, Artificial Chordae. Seminars in Thoracic and Cardiovascular Surgery, Vol 16, No 2 (Summer), 2004: pp 161-168.
- Frater, RWM, Gabbay, S, Shore, D, et al. Reproducible Replacement of Elongated or Ruptured Mitral Valve Chordae. Ann Thorac Surg, 1983; 35; 14-28.
- Gillinov, M, Banbury, M, Pre-Measured Artificial Chordae for Mitral Valve Repair. Ann Thorac Surg, 2007; 84:2127-9.
- Zussa, Claudio, Artificial Chordae in Mitral Valve Surgery. Austin: R.G. Lades Company, 1994: 159.
- Zussa, Claudio, Different Applications of ePTFE Valve Chordae: Surgical Techniques. J Heart Valve Disease, 1996; 5:356-361.

7. ADVERSE REACTIONS

Potential adverse effects associated with the use of any suture include: wound dehiscence, infection, and localized transitory inflammatory tissue reaction.

	2-0 Suture* *Chord-X™ Pre-Measured Loops and Chord-X™ Suture System USP 2-0 suture size is a similar product to Gore-Tex® Suture CV-4. Gore-Tex® is a registered trademark of W. L. Gore & Associates, Inc., Flagstaff, AZ 86004.			
	22mm Taper Point Needle		18mm Taper Point Needle	
Loop Length	3/8 Circle	1/2 Circle	3/8 Circle	1/2 Circle
Adjustable	CXL-20-2238-0	CXL-20-2212-0	CXL-20-1838-0	CXL-20-1812-0
12mm	CXL-20-2238-12	CXL-20-2212-12	CXL-20-1838-12	CXL-20-1812-12
16mm	CXL-20-2238-16	CXL-20-2212-16	CXL-20-1838-16	CXL-20-1812-16
20mm	CXL-20-2238-20	CXL-20-2212-20	CXL-20-1838-20	CXL-20-1812-20
24mm	CXL-20-2238-24	CXL-20-2212-24	CXL-20-1838-24	CXL-20-1812-24

	3-0 Suture* *Chord-X™ Pre-Measured Loops and Chord-X™ Suture System USP 3-0 suture size is a similar product to Gore-Tex® suture CV-5. Gore-Tex® is a registered trademark of W. L. Gore & Associates, Inc., Flagstaff, AZ 86004.			
	22mm Taper Point Needle		18mm Taper Point Needle	
Loop Length	3/8 Circle	1/2 Circle	3/8 Circle	1/2 Circle
Adjustable	CXL-30-2238-0	CXL-30-2212-0	CXL-30-1838-0	CXL-30-1812-0
12mm	CXL-30-2238-12	CXL-30-2212-12	CXL-30-1838-12	CXL-30-1812-12
16mm	CXL-30-2238-16	CXL-30-2212-16	CXL-30-1838-16	CXL-30-1812-16
20mm	CXL-30-2238-20	CXL-30-2212-20	CXL-30-1838-20	CXL-30-1812-20
24mm	CXL-30-2238-24	CXL-30-2212-24	CXL-30-1838-24	CXL-30-1812-24

1 box = 5 sterile single-packaged cases

8. HOW SUPPLIED

All configurations are supplied with (2) 0.118 X 0.276 X 0.073 inches (3 X 7 X 1.85 mm) pledgets. All suture strands used in the construction of the prosthesis are 32 inches in length.

Chord-X™ Pre-Measured Loops and Chord-X™ Suture System for Mitral Chordal Replacement are supplied STERILE unless the integrity of the package has been compromised. Do not use the device if the package has been compromised. The device is for single use only.

9. DEFINITIONS

Use By

Consult Instructions For Use

Do Not Reuse

Catalogue Number

Batch Code

Authorized Representative in the
European Community

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- David, TE, Artificial Chordae.Seminars in Thoracic and Cardiovascular Surgery, Vol 16, No 2 (Summer), 2004: pp 161-168.
- Frater, RWM, Gabbay, S, Shore, D, et al.Reproducible Replacement of Elongated or Ruptured Mitral Valve Chordae.Ann Thorac Surg, 1983; 35: 14-28.
- Gillinov, M, Banbury, M, Pre-Measured Artificial Chordae for Mitral Valve Repair.Ann Thorac Surg, 2007; 84:2127-9.
- Zussa, Claudio, Artificial Chordae in Mitral Valve Surgery.Austin: R.G.Lades Company, 1994: 159.
- Zussa, Claudio, Different Applications of ePTFE Valve Chordae: Surgical Techniques.J Heart Valve Disease, 1996; 5:356-361.

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
FIGURES

Figure 1:

1. Chord-X™ Pre-Measured Loops prosthesis version.

2. Close-up illustration of Chord-X™ Pre-Measured Loops prosthesis version with attached pledget and pre-measured and pre-tied loops.



Figure2:

1. Chord-X™ Suture System prosthesis version.

2. Close-up illustration of Chord-X™ Suture System prosthesis version with attached pledget and untied suture strands.



Figure 3:

1. Removing undesired loop from the Chord-X™ Pre-Measured Loops prosthesis version: Remove undesired prosthetic loops by cutting a **single** strand of suture between the knot and the papillary suture strand.

2. Discard the undesired loop by removing the single strand in entirety.

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

Figure 4:

1. Removing undesired strands from the Chord-X™ Suture System Adjustable prosthesis version: Remove undesired suture strand by cutting a single chordal suture strand.

2. Discard the undesired strand by removing the strand in entirety.

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
Figure 5: Papillary suture placement.


Figure 6: Chord-X™ Pre-Measured Loops leaflet placement and knot tying.


Figure 7: Chord-X™ Suture System Adjustable prosthesis version leaflet placement and knot tying.



Figure 8: Final Chordal Replacement.

On-X Life Technologies, Inc.

1300 East Anderson Lane Building B Austin, Texas 78752 U.S.A.

Telephone: (512) 339-8000 Facsimile: (512) 339-3636

www.onxlti.com • onx@onxlti.com

CryoLife Europa, Ltd.

Bramley House, The Guildway, Old Portsmouth Road, Guildford

Surrey GU3 1LR, United Kingdom

104122 05 Rev F
(07/2017)