## Instructions for Use

### On-X® Ascending Aortic Prosthesis


with the Vascutek Gelweave Valsalva™ Graft

1. English .....	2
2. Bulgaria (Български) .....	18
3. Croatian (Hrvatski) .....	33
4. Czech (Česky) .....	48
5. Danish (Dansk).....	63
6. Dutch (Nederlands).....	78
7. Estonian (Eesti) .....	92
8. Finnish (Suomi).....	107
9. French (Français).....	121
10. German (Deutsch).....	136
11. Greek (Ελληνικά) .....	151
12. Hungarian (Magyar) .....	165
13. Italian (Italiano).....	180
14. Japanese (日本語).....	194
15. Korean (한국어).....	209
16. Latvian (Latviski) .....	223
17. Lithuanian (Lietuvių k.) .....	237
18. Norwegian (Norsk).....	252
19. Persian (پارسی).....	267
20. Polish (Polski).....	282
21. Portuguese (Português).....	297
22. Romanian (Română).....	312
23. Russian (На русском языке) .....	326
24. Serbian (Српски) .....	340
25. Slovak (Slovensky) .....	355
26. Slovenian (Slovenščina) .....	370
27. Spanish (Español) .....	385
28. Swedish (Svenska) .....	400
29. Turkish (Türkçe).....	415
30. Ukrainian (Українська).....	429
Figures and Tables.....	442

# On-X® Ascending Aortic Prosthesis

with the Vascutek Gelweave Valsalva™ Graft

## Instructions for Use

Ascending Aortic Prosthesis

On-X® Ascending Aortic Prosthesis REF ONXAAP

# TABLE OF CONTENTS

<b>1. DEVICE DESCRIPTION</b> .....	<b>4</b>
<b>2. INDICATIONS</b> .....	<b>4</b>
<b>3. CONTRAINDICATIONS</b> .....	<b>4</b>
<b>4. WARNINGS AND PRECAUTIONS</b> .....	<b>5</b>
4.1 Warnings.....	5
4.2 Precautions .....	6
<b>5. POTENTIAL ADVERSE EVENTS</b> .....	<b>7</b>
<b>6. INDIVIDUALIZATION OF TREATMENT</b> .....	<b>7</b>
6.1 Specific Patient Population.....	8
<b>7. PATIENT COUNSELING</b> .....	<b>8</b>
<b>8. HOW SUPPLIED</b> .....	<b>8</b>
8.1 Available Models and Sizes .....	8
8.2 Ascending Aortic Prosthesis Holder Handle .....	8
8.3 Packaging .....	8
8.4 Storage .....	9
8.5 Accessories .....	9
<b>9. DIRECTIONS FOR USE</b> .....	<b>10</b>
9.1 Physician Training .....	10
9.2 Sterilization and Re-sterilization.....	10
9.3 Handling and Preparation Instructions .....	10
9.4 Device Implantation.....	12
9.5 Valve Suturing Techniques.....	13
9.6 Graft Length and Suturing.....	13
9.7 Leaflet Motion Assessment and Valve Rotation .....	14
9.8 Valve Orientation .....	14
<b>10. POSTOPERATIVE INFORMATION</b> .....	<b>15</b>
10.1 Magnetic Resonance Imaging (MRI) Compatibility .....	15
10.2 Returned Goods .....	16
<b>11. PATIENT INFORMATION</b> .....	<b>16</b>
11.1 Patient Registration .....	16
11.2 Patient Record Card.....	16
<b>12. DISCLAIMER OF WARRANTIES</b> .....	<b>16</b>

**FIGURES AND TABLES..... 442**

Figure 1. On-X® Ascending Aortic Prosthesis .....442

Figure 2. On-X® Ascending Aortic Prosthesis with  
Holder Handle Assembly.....442

Figure 3. Replica Sizer .....443

Figure 4. Instrument Handle .....443

Figure 5. Valve Rotator .....443

Figure 6. Leaflet Probe .....443

Figure 7. Removing Prosthesis from Inner Package .....443

Figure 8. Removing Holder Handle from Prosthesis .....444

Table 1. On-X Ascending Aortic Prosthesis Specifications (mm).....445

Table 2. Definitions .....446

**Please note: All Figures and Tables are displayed in the last pages of this booklet. Please refer to the final pages for illustrations and data.**

## **1. DEVICE DESCRIPTION**

Manufactured by On-X Life Technologies, Inc.<sup>TM</sup> (On-X LTI), the On-X<sup>®</sup> Ascending Aortic Prosthesis combines the On-X<sup>®</sup> Prosthetic Heart Valve and the Gelweave Valsalva<sup>TM</sup> Vascular Prosthesis (Figure 1). The prosthesis is available in sizes 19, 21, 23, 25, and 27/29 mm. The Gelweave Valsalva<sup>TM</sup> Vascular Prosthesis is incorporated into the cuff structure of the On-X valve to create the ascending aortic prosthesis. Orientation reference marks are provided on the exterior of the graft to aid implantation.

The On-X<sup>®</sup> Prosthetic Heart Valve is a bileaflet mechanical heart valve, which consists of an orifice and two leaflets. The orifice inflow area has a flared inlet designed to reduce flow turbulence, and the outflow rim consists of leaflet guards designed to protect the leaflets while in the closed position. The leaflets rotate around tabs located within the inner circumference of the orifice ring. In the closed position, each leaflet forms a nominal angle of 40° relative to the plane of the orifice. In the open position, the plane of each leaflet forms a nominal angle of 90° relative to the plane of the orifice. The leaflets have a travel arc of 50° to the closed position. The orifice is composed of a graphite substrate coated with On-X<sup>®</sup> Carbon, a pure unalloyed form of pyrolytic carbon. The leaflets consist of On-X<sup>®</sup> Carbon deposited on graphite substrates, which are impregnated with 10 weight% tungsten to provide radiopacity. The sewing cuff is constructed of polytetrafluoroethylene (PTFE) fabric mounted on the orifice using titanium retaining rings and 5-0 polyester suture material. This form of sewing cuff attachment to the orifice allows for rotation of the sewing cuff in situ during implantation.

The Gelweave Valsalva<sup>TM</sup> Vascular Prosthesis is a woven polyester prosthesis designed to mimic the geometry of the sinus of Valsalva. The graft is impregnated with an absorbable protein (modified mammalian gelatin) to eliminate preclotting of the graft in surgery. The gelatin has been cross-linked to a set level to control its rate of removal. After implant, the gelatin is hydrolyzed within approximately 14 days and is replaced by normal tissue incorporation. Gelatin is a non-antigenic and non-toxic protein as reflected in its use as a safe plasma expander.

The gelatin used in the graft is manufactured from animals native to and exclusively raised in Australia.

## **2. INDICATIONS**

The On-X<sup>®</sup> Ascending Aortic Prosthesis is indicated for the replacement of a diseased, damaged, or malfunctioning native or prosthetic heart valve in the aortic position in cases that involve an ascending aortic aneurysm or other associated aortic disease.

## **3. CONTRAINDICATIONS**

The On-X<sup>®</sup> Ascending Aortic Prosthesis is contraindicated for patients unable to tolerate anticoagulation therapy or who exhibit sensitivity to polyester or materials of bovine origin.

## **4. WARNINGS AND PRECAUTIONS**

### **4.1 Warnings**

**USE ONLY AS INDICATED** (see section 2).

#### **FOR SINGLE USE ONLY**

**DO NOT** use the On-X<sup>®</sup> Ascending Aortic Prosthesis if:

- the prosthesis has been dropped, damaged, or mishandled in any way;
- the expiration date has elapsed;
- the tamper evident seals are broken;
- package damage is evident;
- the serial number tag does not match the serial number on the container label.

**DO NOT** pass a catheter, surgical instrument, or transvenous pacing lead through the aortic prosthetic heart valve as this may cause valvular insufficiency, leaflet damage, leaflet dislodgment, and/or catheter/ instrument/ lead entrapment.

**DO NOT** re-sterilize the On-X<sup>®</sup> Ascending Aortic Prosthesis.

**DO NOT** use the On-X<sup>®</sup> Ascending Aortic Prosthesis if the expiration date has elapsed. If a prosthesis is unused, its foil pouch and sealed plastic tray container is undamaged, and the sterility expiration date has passed, the prosthesis should be returned to On-X LTI.

**DO NOT** use the On-X<sup>®</sup> Ascending Aortic Prosthesis if either tamper evident label at each end of the outer box has been broken or is damaged. If either tamper evident label has been broken or is damaged, use another prosthesis and arrange to return the unused product to On-X LTI.

**DO NOT** use the On-X<sup>®</sup> Ascending Aortic Prosthesis if the protective foil pouch is damaged or if the prosthesis has been dropped, damaged, or mishandled in any way. If any damage is found, use another prosthesis and arrange to return the unused product to On-X LTI.

**DO NOT** preclot the Gelweave Valsalva<sup>™</sup> Vascular Prosthesis.

**DO NOT** exert excessive force on the valve components or the graft.

**DO NOT** attempt to re-insert the holder handle into the prosthesis (valve) once it has been released from the prosthesis. Doing so could damage the valve and cause catastrophic failure of the valve or harm to the patient.

**DO NOT** use the holder handle to rotate the valve once implanted.

**DO NOT** attempt to disassemble the holder handle. A small spring protected inside the assembled holder handle could escape and be unnoticed in the sterile field.

**ONCE** the foil package is opened, the prosthesis must be implanted within 24 hours. If the plastic trays with sealed Tyvek<sup>®</sup> lids are not damaged, sterility of the prosthesis is maintained and the sealed tray set can remain out of the foil pouch for up to 24 hours. If the prosthesis is not implanted within 24 hours of opening the foil pouch, arrange return of the prosthesis to On-X LTI.

**TEST** the leaflet mobility only with the On-X leaflet probe.

All accessory instruments must be cleaned and sterilized prior to use according to the instructions provided separately with the instruments.

**USE ONLY** the On-X rotator to rotate the valve in situ. Use only the correspondingly sized rotator. Use of the wrong size rotator could damage the valve.

**DO NOT** use the On-X® Ascending Aortic Prosthesis if the inner tray lid is damaged. If the inner tray lid is damaged, use another prosthesis and arrange to return the unused product to On-X LTI

**DO NOT** use the On-X® Ascending Aortic Prosthesis if the valve does not rotate easily. Use another prosthesis and arrange to return the unused product to On-X LTI.

**DO NOT** use the On-X® Ascending Aortic Prosthesis if the serial number tag does not match the container label. Use another prosthesis and arrange to return the unused product to On-X LTI.

**DO NOT** size the sewing cuff of the size 19 mm through 25 mm aortic valve to fit within the annulus. Size 27/29 mm aortic valves are designed to be placed in an intra-annular position and have a replica sizer to mimic this placement.

#### 4.2 Precautions

**CAUTION:** Federal (USA) law restricts this device to sale, distribution and use by or on the order of a physician.

Handle the prosthesis with On-X Valve Instruments and the provided product handle. Only On-X valve sizers should be used during the selection of the prosthesis size; other sizers may result in improper prosthesis selection.

Avoid contacting the carbon surfaces of the valve with gloved fingers or any metallic or abrasive instruments as they may cause damage to the valve surface not seen with the unaided eye that may lead to accelerated valve structural dysfunction, leaflet escape, or serve as a nidus for thrombus formation.

No resistance should be experienced when inserting the rotator. If resistance is encountered, stop, remove, and realign the rotator before attempting to insert the rotator again. Retest leaflet motion after rotation. If free leaflet motion cannot be achieved, remove the prosthesis.

Clamping may damage any vascular graft. Handling of the graft must be gentle and minimized in order to prevent loss of gelatin coating.

Excessive tension on the graft should be avoided.

The use of round taper point needles minimizes graft damage. When de-airing through the graft, a 19-gauge needle is usually sufficient. Cutting point needles should be avoided to prevent damage.

Gelweave Valsalva™ Vascular Prosthesis has a woven structure and therefore should be cut with a cautery to minimize fraying. Use of the cautery for any sealed graft can cause burning. Immersion of the graft in saline prior to use of the cautery will prevent focal burning. Immersion of the ascending aortic prosthesis should be no longer than 5 minutes. The prosthesis should not be allowed to dry out after soaking.

Do not attempt to rotate the valve if any significant resistance to rotation is encountered and only use the valve rotator attached to the instrument handle to rotate the valve. The torque required to rotate the valve in situ should be about the same as that required when testing rotation before implantation. If noticeably greater torque is required to rotate, stop attempting rotation. If rotation is necessary and cannot be performed, remove the prosthesis.

Use only On-X valve sizers when sizing the annulus. Sizers contain cylindrical, conical, and aortic replica ends.

## **5. POTENTIAL ADVERSE EVENTS**

Adverse events potentially associated with the use of aortic valved grafts include, but are not limited to:

- angina
- cardiac arrhythmia
- endocarditis
- heart failure
- hemolysis
- hemolytic anemia
- hemorrhage
- myocardial infarction
- valve leaflet entrapment (impingement)
- valve nonstructural dysfunction
- valve pannus overgrowth
- valve perivalvular leak
- valve regurgitation
- valve structural dysfunction
- valve thrombosis
- thromboembolism
- stroke
- graft infection
- false aneurysm, aneurysm
- secondary fistula formation
- graft occlusion or kinking
- excessive blood loss through the graft

It is possible that these complications could lead to:

- reoperation
- permanent disability
- explantation
- death

Mechanical prosthetic heart valves produce audible sounds as a normal function of their operation. In some patients, these sounds may be objectionable.

Risk of Re-use:

In accordance with the EU Medical Device Directive, 93/42/EEC, Annex I, Section 13.6h, the device manufacturer must provide information on risks associated with re-use of a single use device. Therefore, the following statement is provided:

Risks of re-using the AAP device:

The implanted graft is customized to the patient. In addition to the risks listed in this section, re-use in another patient, would risk exsanguination. Additional risks are bloodborne pathogens, and bacterial infection.

## **6. INDIVIDUALIZATION OF TREATMENT**

Adequate anticoagulant or anticoagulant/antiplatelet therapy should be administered postoperatively. Selection of an anticoagulant or anticoagulant/antiplatelet regimen is based on the particular needs of the patient, the clinical situation and established professional guidelines.

## 6.1 Specific Patient Population

The safety and effectiveness of the On-X<sup>®</sup> Ascending Aortic Prosthesis has not been established for the following specific populations because it has not been studied in these populations:

- patients who are pregnant;
- nursing mothers;
- patients with chronic endocarditis;
- patients requiring pulmonary root replacement.

## 7. PATIENT COUNSELING

- Prophylactic antibiotic treatment must be provided to all patients with prosthetic valves undergoing dental procedures or other potentially bacteremic procedures.
- Patients require anticoagulation or anticoagulant/antiplatelet therapy.
- Patients should be encouraged to complete the Patient Record Card provided with the prosthesis and carry it with them at all times.
- Patients should be informed of the use of bovine-derived materials within this device in advance of procedure.

## 8. HOW SUPPLIED

### 8.1 Available Models and Sizes

The On-X<sup>®</sup> Ascending Aortic Prosthesis is available in sizes 19, 21, 23, 25, and 27/29 mm. The symbol SZ mm on the box and container labels refers to the tissue annulus diameter of the valve in millimeters. A circular symbol representing actual graft diameter appears on each box end label. See Table 1 for all product specifications.

### 8.2 Ascending Aortic Prosthesis Holder Handle

The On-X<sup>®</sup> Ascending Aortic Prosthesis is supplied mounted on a plastic extended holder handle with spring action for one time use. This holder handle (Figure 2) facilitates positioning the prosthesis during surgery. The prosthesis can be released from the holder by depressing the central plunger.

**WARNING:** Do not attempt to re-insert the holder handle into the prosthesis (valve) once it has been released from the prosthesis. Doing so could damage the valve and cause catastrophic failure of the valve or harm to the patient.

**WARNING:** Do not attempt to disassemble the holder handle. A small spring protected inside the assembled holder handle could escape and be unnoticed in the sterile field.

### 8.3 Packaging

The On-X<sup>®</sup> Ascending Aortic Prosthesis is provided sterile, mounted on a holder, in inner and outer plastic trays enclosed in a foil pouch that serves as a vapor barrier and preserves optimal prosthesis characteristics. The foil

pouch and outer tray are not sterile. Tyvek® seals on the inner and outer trays must be intact. Any damage to the trays renders the prosthesis non-sterile. In the event of damage to the primary packaging, the product must not be used and should be returned to On-X LTI. The packaging consists of the following items:

- Outer box with labeling
- Plastic prosthesis tray x 2, inner and outer
- Plastic extended valve holder with handle
- Prosthesis serial number tag
- Silica gel sachet
- Foil pouch
- May contain Instructions for Use booklet
- Patient record card
- Instructional card, complete with Courtesy Labels, used for accessing: 1) Electronic Implant Registration and 2) Electronic Instructions for Use

Instruments for implantation of the On-X® Ascending Aortic Prosthesis are supplied separately, NON-STERILE, and must be cleaned and sterilized prior to use as outlined in the Instructions for Use provided with the instruments.

#### 8.4 Storage

The On-X® Ascending Aortic Prosthesis is provided sterile. The sterility expiration date of the device is recorded on the outer package labeling. Appropriate inventory control should be maintained so that prostheses with earlier expiration dates are preferentially implanted and expiration is avoided. To protect the prosthesis, it should be stored in its outer box until used. The storage environment should be clean, cool, and dry at a temperature not less than 0°C (32° F) or higher than 35°C (95° F).

**WARNING:** Once the foil package containing the sealed plastic tray set and prosthesis is opened, the prosthesis must be implanted within 24 hours. If the plastic trays with sealed Tyvek® lids are not damaged, sterility of the prosthesis is maintained and the sealed tray set can remain out of the foil pouch for up to 24 hours. If the prosthesis is not used within 24 hours, contact On-X LTI Customer Service (+1 512-339-8000 or cs@onxlti.com) to arrange return of the prosthesis.

#### 8.5 Accessories

The On-X® Ascending Aortic Prosthesis (valve) is designed to be used only with On-X instruments. The instruments, supplied separately, are provided as a set, which includes sizers, rotators, an instrument handle, and a leaflet probe. The instruments are reusable. See the On-X Prosthetic Heart Valve Instruments Instructions for Use for complete descriptions of instruments, their use and reprocessing.

##### **Replica Sizers**

Aortic replica sizers are provided for all aortic valve sizes (Figure 3). They are used to assure fit of the aortic valve without obstruction of the coronary arteries. The size 19 through 25 aortic replica sizer shape is intended to model supra-annular positioning of the valve sewing ring. Size 27/29 aortic replica sizer is intended to model intra-annular positioning.

##### **Instrument Handle**

The instrument handle is comprised of a grip, a bendable stem and a tip that facilitates holding the rotator during surgery (Figure 4).

**Rotator**

The valve rotator (Figure 5) is used for reorienting an in situ valve after the conduit is cut to length. The rotator consists of a plastic head with a centrally located leaflet probe and a slot for connection to the instrument handle. To attach the rotator to the instrument handle, insert the instrument handle tip directly into the slot on the end of the rotator handle. The rotator snaps into place after the application of a light insertion force.

**Leaflet Probe**

The leaflet probe (Figure 6) is a flexible rod with tapered ends. The leaflet probe may be used to gently move the leaflets to verify that they open and close freely.

**9. DIRECTIONS FOR USE**

**WARNING:** Do not use the On-X® Ascending Aortic Prosthesis if:

- the prosthesis has been dropped, damaged, or mishandled in any way;
- the expiration date has elapsed;
- the tamper evident seal is broken;
- package damage is evident;
- the serial number tag does not match the serial number on the container label.

**9.1 Physician Training**

The techniques for implanting the On-X® Ascending Aortic Prosthesis are similar to those used for any ascending aortic prosthesis. Trained cardiovascular surgeons will be familiar with these techniques.

**9.2 Sterilization and Re-sterilization**

The On-X® Ascending Aortic Prosthesis is provided sterile. If the sterility expiration date has passed, or if upon removal from the outer box, the foil pouch is punctured or the sterility barriers of the sealed plastic trays inside the foil pouch are compromised, do not use the prosthesis, contact On-X LTI Customer Service (+1 512-339-8000 or [cs@onxlti.com](mailto:cs@onxlti.com)) to arrange return of the prosthesis.

**WARNING:** Once the foil package is opened, the prosthesis must be implanted within 24 hours. If the plastic trays with sealed Tyvek® lids are not damaged, sterility of the prosthesis is maintained and the sealed tray set can remain out of the foil pouch for up to 24 hours. If the prosthesis is not implanted within 24 hours of opening the foil pouch, contact On-X LTI Customer Service (+1 512-339-8000 or [cs@onxlti.com](mailto:cs@onxlti.com)) to arrange return of the prosthesis.

**WARNING:** Do not re-sterilize the On-X® Ascending Aortic Prosthesis.

**9.3 Handling and Preparation Instructions**

**CAUTION:** Handle the prosthesis only with On-X LTI Instruments. Only On-X valve sizers should be used during the selection of the prosthesis size; other sizers may result in improper prosthesis selection.

**CAUTION:** Avoid contacting the carbon surfaces of the valve with gloved fingers or any metallic or abrasive instruments as they may cause damage to the valve surface not seen with the unaided eye that may lead to accelerated valve structural dysfunction, leaflet escape, or serve as a nidus for thrombus formation.

**CAUTION:** Avoid damaging the prosthesis through the application of excessive force to the valve orifice or leaflets or to the graft.

### **Circulating Nurse**

1. Check the expiration date on the outer box.

**WARNING:** Do not use the On-X® Ascending Aortic Prosthesis if the expiration date has elapsed. If a prosthesis is unused, its foil pouch and sealed plastic tray container is undamaged, and the sterility expiration date has passed, the prosthesis should be returned to On-X LTI.

2. Check the labels at each end of the box for integrity.

**WARNING:** Do not use the On-X® Ascending Aortic Prosthesis if either tamper evident label at each end of the outer box has been broken or is damaged. If either tamper evident label has been broken or is damaged, use another prosthesis and arrange for return of the unused product through On-X LTI Customer Service (+1 512-339-8000 or [cs@onxlti.com](mailto:cs@onxlti.com)).

3. Open the outer box and remove the protective foil pouch containing the prosthesis and package inserts. Inspect the foil pouch for damage.

**WARNING:** Do not use the On-X® Ascending Aortic Prosthesis if the protective foil pouch is damaged or if the prosthesis has been dropped, damaged, or mishandled in any way. If any damage is found, use another prosthesis and arrange for return of the unused product through On-X LTI Customer Service (+1 512-339-8000 or [cs@onxlti.com](mailto:cs@onxlti.com)).

4. Complete electronic implant registration on [www.onxlti.com/implantregistration](http://www.onxlti.com/implantregistration) as completely as local law allows and as soon as possible. This allows the patient to be entered into the tracking database, which could be important for future notices regarding the prosthesis. Give the patient record card to the patient.
5. Open the foil pouch by tearing it at one of the notches on either side of the indicated tear line. Illustrations for opening the foil pouch, removing and opening the outer tray, introducing the inner tray into the sterile field and opening the inner tray are given on the foil pouch.
6. Open the outer tray lid by grasping the pull tab marked on the corner of the outer lid and pulling the lid in the direction of the arrow. (Refer to the illustration on the foil pouch.)
7. The inner container can be placed in the sterile field by gently inverting the outer container slightly above the sterile field and allowing the inner container to slip out into the sterile field. (Refer to the illustration on the foil pouch.)

**Scrub Nurse/Surgeon:**

8. The scrub nurse may open the sterile inner tray by grasping the pull tab marked on the corner of the inner lid and pulling toward the direction of the arrow. (Refer to the illustration on the foil pouch.) The prosthesis inner container is then placed onto the instrument tray.

**WARNING:** Do not use the On-X® Ascending Aortic Prosthesis if the inner tray lid is damaged. If the inner tray lid is damaged, use another prosthesis and arrange for return of the unused product through On-X LTI Customer Service (+1 512-339-8000 or cs@onxlti.com).

9. To remove the prosthesis from the inner tray, stabilize the tray on the sterile field with one hand. With the other hand, insert the fingers around the holder handle end and gently lift upward with enough force to release it from the tray (Figure 7, and the illustration on the inner tray lid).
10. To test rotation, carefully grasp the sewing cuff with a gloved hand using a light grip and gently turn the top of the holder handle in either direction. The valve should easily rotate within the sewing cuff.

**WARNING:** Do not use the On-X® Ascending Aortic Prosthesis if the valve does not rotate easily. Use another prosthesis and arrange for return of the unused product through On-X LTI Customer Service (+1 512-339-8000 or cs@onxlti.com).

11. Check the serial number tag against the label on the outer tray.

**WARNING:** Do not use the On-X® Ascending Aortic Prosthesis if the serial number tag does not match the container label. Use another prosthesis and arrange for return of the unused product through On-X LTI Customer Service (+1 512-339-8000 or cs@onxlti.com).

12. Remove the serial number tag by cutting the suture that holds it on the prosthesis. If desired, the tag can be used to check for sterility by standard culture techniques immediately after it is removed. The prosthesis is now ready for implant.

**9.4 Device Implantation**

**WARNING:** All accessory instruments must be cleaned and sterilized prior to use according to the On-X Instrumentation Instructions for Use provided separately with the instruments.

**Sizing**

**CAUTION:** Use only On-X valve sizers when sizing the annulus. Sizers contain cylindrical, conical, and aortic replica ends.

Cylindrical sizers correspond to valve sizes 19 mm through 25 mm. Conical sizers correspond to valve size 27/29 mm. The correct valve size is determined by obtaining a comfortable, not tight, fit of the sizer within the annulus. When a comfortable fit is found, the corresponding valve size is signified by the identification on the sizer.

Aortic replica sizers are provided for all aortic valve sizes. For size 19 mm through 25 mm aortic valves, the aortic replica sizers are used to verify that the aortic valve can be properly seated in the annulus and that the coronary arteries remain unobstructed. Size 19 mm through 25 mm aortic valves are designed to fit within the annulus at implant such that the exposed carbon flare rests in the annulus and the sewing ring is intra-supra annular.

**WARNING:** Do not size the sewing cuff of the size 19 mm through 25 mm aortic valve to fit within the annulus. Size 27/29 mm aortic valves are designed to be placed in an intra-annular position and have a replica sizer to mimic this placement.

## 9.5 Valve Suturing Techniques

Valve suturing techniques vary according to the preferences of the implanting surgeon and patient condition. The aortic valve is designed to have the tissue annulus about the orifice flare. The general consensus among surgeons is that a non-everting interrupted mattress suture technique, with or without pledgets, provides the best conformation of the valve annulus to the outer surface of the flare.

Sutures should be passed through the midpoint of the sewing cuff. This allows the sewing cuff to remain flexible and conform to the annulus. It also prevents the suture needle from contacting the titanium rings that lie within the sewing cuff. The orientation marks on the sewing cuff may be used to aid in suture placement.

When all sutures are in place, the prosthesis is advanced into the annulus and the sutures are tied down. It is suggested that the first 3 knots be tied equidistant to one another and midway between the commissures to stabilize the valve in the annulus. The holder handle is removed from the valve by carefully depressing the plunger in the holder handle as shown in Figure 8, then gently lifting the holder handle out of the valve/prosthesis.

**WARNING:** Do not attempt to reinsert the valve holder handle into the valve once it has been removed. Doing so could damage the valve and cause catastrophic failure of the valve or harm to the patient.

## 9.6 Graft Length and Suturing

The graft portion of the prosthesis may be trimmed at any time during the implant procedure to the appropriate length.

**CAUTION:** Clamping may damage any vascular graft. Handling of the graft must be gentle and minimized in order to prevent loss of gelatin coating.

**CAUTION:** Excessive tension on the graft should be avoided.

**CAUTION:** The use of round taper point needles minimizes graft damage. When de-airing through the graft, a 19-gauge needle is usually sufficient. Cutting point needles should be avoided to prevent damage.

**CAUTION:** Gelweave Valsalva™ Vascular Prosthesis is a woven structure and therefore should be cut with a cautery to minimize fraying. Immersion of the graft portion of the prosthesis in saline immediately prior to use will prevent focal burning which may result during cauterization for coronary ostia or to trim the graft. The prosthesis should be immersed in saline for no longer than 5 minutes and should not be allowed to dry out after soaking.

## 9.7 Leaflet Motion Assessment and Valve Rotation

### **Leaflet Motion Testing**

Once the prosthesis is in place, free motion of the leaflets must be tested. To test leaflet mobility, use the leaflet probe to gently move the leaflets to verify that they open and close freely.

**WARNING:** Test the leaflet mobility only with the On-X LTI leaflet probe.

### **Rotation**

If the leaflets do not move freely, gently rotate the valve in either direction until it reaches a position where leaflet interference is not encountered.

The rotator may be used with or without the instrument handle attached. As needed, attach the instrument handle to the rotator by inserting the instrument handle tip into the slot on the end of the rotator handle until it snaps firmly into position. With the rotator leaflet probe between the leaflets and the cross-bar aligned with the leaflet pivot axis of the valve, carefully insert the valve rotator into the valve until it seats easily in place.

**CAUTION:** Do not attempt to rotate the valve if any significant resistance to rotation is encountered. The torque required to rotate the valve in situ should be about the same as that required when testing rotation before implantation. If noticeably greater torque is required to rotate, stop attempting rotation. If rotation is necessary and cannot be performed, remove the prosthesis.

**WARNING:** Do not use the holder handle to rotate the valve once implanted.

**WARNING:** Use only the On-X rotator to rotate the valve in situ. Use only the correspondingly sized rotator. Use of the wrong size rotator could damage the valve.

**CAUTION:** No resistance should be experienced when inserting the rotator. If resistance is encountered, stop, remove, and realign the rotator before attempting to insert the rotator again. Retest leaflet motion after rotation. If free leaflet motion cannot be achieved, remove the prosthesis.

## 9.8 Valve Orientation

Based on clinical studies, there is no preferred orientation for the On-X® Prosthetic Heart Valve in the On-X® Ascending Aortic Prosthesis.

## 10. POSTOPERATIVE INFORMATION

### 10.1 Magnetic Resonance Imaging (MRI) Compatibility

#### MR Conditional

The **On-X Prosthetic Heart Valve, Conform-X Mitral Heart Valve Prosthesis, Size 25-33\***, was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-08. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania.

Non-clinical testing demonstrated that the On-X Conform-X Mitral Heart Valve Prosthesis, Size 25-33, is MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

#### Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

#### MRI-Related Heating

In non-clinical testing, the On-X Prosthetic Heart Valve, Conform-X Mitral Heart Valve Prosthesis, Size 25-33, produced the following temperature rise during MRI performed for 15-min of scanning (i.e., per pulse sequence) in the 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR system:

*Highest temperature change +1.6°C*

Therefore, the MRI-related heating experiments for the On-X Conform-X Mitral Heart Valve Prosthesis, Size 25-33, at 3-Tesla, using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9-W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7-W/kg), indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.6°C.

#### Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the On-X Conform-X Mitral Heart Valve Prosthesis, Size 25-33. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	1,090-mm <sup>2</sup>	686-mm <sup>2</sup>	1,478-mm <sup>2</sup>	1,014-mm <sup>2</sup>
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

\*The MRI findings apply to this specific largest heart valve prosthesis and all other smaller sizes which are made from similar materials.

## **10.2 Returned Goods**

Prior authorization from On-X LTI Customer Service is required for the return of any product. For any questions regarding the valve or for return authorization, please contact Customer Service at +1 512-339-8000 or [cs@onxlti.com](mailto:cs@onxlti.com).

Licensed under U.S. Patent Nos. 5,137,532; 5,545,216; 5,772,694; 5,908,452; 5,677,061; 6,096,075; 5,641,324; and foreign counterpart patents.

## **11. PATIENT INFORMATION**

### **11.1 Patient Registration**

In each On-X® Ascending Aortic Prosthesis package, an Instructional Card, complete with Courtesy Labels, used for accessing the Electronic Implant Registration and Electronic Instructions for Use is included. On-X LTI requests that the Electronic Implant Registration be completed immediately.

On-X LTI will use these data for notification purposes and to help with inventory restocking in the hospital. All patient information remains strictly confidential, and the release of patient identifying information can be refused if allowed by local law.

### **11.2 Patient Record Card**

A Patient Record Card is provided with the prosthesis. The card should be given to the patient (or placed in the patient record to be given to the patient). Patients should be encouraged to complete the card and carry it with them at all times.

## **12. DISCLAIMER OF WARRANTIES**

Because of the complications listed previously that may occur with the use of any ascending aortic prosthesis and the possibilities of damage, also noted previously, before, during or after implantation, On-X LTI warrants only that the product shall conform to On-X LTI's standard specifications. No other warranty is made by On-X LTI concerning the function of the product in use, and On-X LTI assumes no risk whatsoever as to the results of the use of this product. The entire risk with use of the product is that of the buyer. On-X LTI disclaims all other warranties, respecting the product, expressed or implied, including but not limited to those related to the product's merchantability or fitness for a particular purpose. On-X LTI shall not be liable for any direct, special, consequential or incidental loss, damage or expense related to the use of the product. No person has any authority to alter any of these conditions or to bind On-X LTI to any additional responsibility or warranty in connection with the use of the product.

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Figure 1. On-X® Ascending Aortic Prosthesis





Figure 2. On-X® Ascending Aortic Prosthesis with Holder Handle Assembly


**Figure 3. Replica Sizer**

**Figure 4. Instrument Handle**

**Figure 5. Valve Rotator**

**Figure 6. Leaflet Probe**

**Figure 7. Removing Prosthesis from Inner Package**



Figure 8. Removing Holder Handle from Prosthesis

**Table 1. On-X Ascending Aortic Prosthesis Specifications (mm)**

Model Designator	Valve Size/Type	Tissue Annulus (mounting) Diameter	Graft Internal Diameter	Valve Orifice Internal Diameter	External Sewing Cuff Diameter	Valve Orifice Internal Area (mm <sup>2</sup> )
ONXAAP-19	19 Aortic	19	20	17.4	27	228
ONXAAP-21	21 Aortic	21	22	19.4	30	284
ONXAAP-23	23 Aortic	23	24	21.4	33	344
ONXAAP-25	25 Aortic	25	26	23.4	34	411
ONXAAP-27/29	27/29 Aortic	27-29	26	23.4	36	411

**Table 2. Definitions**

	Authorized representative in the European community		Do not use if package is damaged
	Manufacturer		Do not re-use
	Consult instructions for use		Use-by date
 <small>www.onxlti.com/ifu/AAP</small>	Consult instructions for use		Serial number
	Catalogue number		Date of manufacture
	Sterilized using Ethylene Oxide		Do not re-sterilize
	MR conditional		Temperature limitation
	Federal Law (USA) restricts this device to sale by or on the order of a physician.		Product must be used within 24 hours of opening foil pouch.
	Graft Diameter in millimeters		Do not discard
	Caution, consult accompanying documents	SZ mm	Valve Size in millimeters
	The On-X Prosthetic Heart Valve was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials International, Designation: F2503-08. See <a href="http://www.onxlti.com/onxlti-hvm-mri-memo-3-0-tesla-march-2010.html">www.onxlti.com/onxlti-hvm-mri-memo-3-0-tesla-march-2010.html</a> for more information.		

On-X Life Technologies, Inc.

1300 East Anderson Lane, Bldg. B, Austin, Texas 78752 U.S.A.

TEL: (512) 339-8000 • FAX: (512) 339-3636

www.onxlti.com • onx@onxlti.com

EC

REP

CryoLife Europa, Ltd.

Bramley House, The Guildway, Old Portsmouth Road

Guildford, Surrey GU3 1LR, United Kingdom