Prosthetic Heart Valve Instruments Instructions for Use

1. English	1
2. Bulgaria (БЪЛГАРСКИ)	9
3. Brazilian Portuguese (PORTUGUESE)	17
4. Czech (ČESKY)	25
5. Danish (DANSK)	33
6. German (DEUTSCH)	41
7. Greek (ΕΛΛΗΝΙΚΑ)	49
8. Spanish (ESPAÑOL)	57
9. Estonian (EESTI)	65
10. Finnish (SUOMI)	73
11. French (FRANÇAIS)	81
12. Croatian (HRVATSKI)	89
13. Hungarian (MAGYAR)	97
14. Indonesia (INDONESIA)	105
15. Italian (ITALIANO)	113
16. Japanese (日本語)	121
17. Korean (한국어)	129
18. Lithuanian (LIETUVIŲ K.)	137
19. Latvian (LATVIEŠU)	145
20. Dutch (NEDERLANDS)	153
21. Norwegian (NORSK)	161
22. Polish (POLSKI)	169
23. Portuguese (PORTUGUÊS)	177
24. Romanian (ROMÂNĂ)	185
25. Russian (РУССКИЙ)	193
26. Slovak (SLOVENČINA)	201
27. Slovenian (SLOVENSKO)	209
28. Serbian (СРПСКИ)	217
29. Swedish (SVENSKA)	225
30. Turkish (TÜRKÇE)	233
31. Ukrainian (УКРАЇНСЬКА)	241
32. Chinese (中文)	249
33. Persian (فارسی)	257

On-X®

Prosthetic Heart Valve Instruments

Инструменти за сърдечна клапна протеза
Instrumentos para prótese valvular cardíaca
Nástroje k protetické srdeční chlopní
Instrumenter til hjerteklapprotese
Instrumente für Herzklappenprothesen
Εργαλεία προθετικής καρδιακής βαλβίδας
Instrumental para Válvula Cardíaca Protésica
Südameklapiproteesi instrumendid
Sydänläppäproteesien asennusvälineet
Instruments de prothèse valvulaire cardiaque
Instrumenti za protetske srčane zaliske
Szívbillentyű protézis műszerek
Instrumen Katup Jantung Prostetik
Strumentario per protesi valvolare cardiaca
人工心臓弁インストルメント
인공심장판막 기구

Protezinio širdies vožtuvo instrumentai
Sirds vārstulja protēzes instrumenti
Instrumenten voor hartklepprothesen
Hjerteklapproteeseinstrumenter
Instrumenty do protezy zastawki serca
Instrumentos para Prótese Valvular Cardíaca
Instrumente pentru proteze valvulare cardiace
Инструменты для протезов клапанов сердца
Nástroje na manipuláciu so srdcovou chlopňovou protézou
Instrumenti za protetično srčno zaklopko
Инструменти за протетске срчане залиске
Instrument för konstgjord hjärtklaff
Prostetik Kalp Kapağı Aletleri
Инструменти для протезів клапана серця
人工心脏瓣膜器械
ایزار دریچه قلب مصنوعی

Instructions for Use

Инструкция за употреба
Instruções de uso
Návod k použití
Bruksanvisning
Gebrauchsanweisung
Οδηγίες χρήσης
Mode de Empleo
Kasutusjuhend
Käyttöohje
Mode d'emploi
Upute za uporabu
Használati utasítás
Instruksi Penggunaan
Istruzioni per l'uso
取扱説明書
사용설명서

Naudojimo instrukcijos
Lietošanas norādījumi
Gebruiksaanwijzing
Bruksanvisning
Instrukcja użytkowania
Instruções de utilização
Instrucţiuni de utilizare
Инструкции по применению
Návod na použitie
Navodila za uporabo
Упутство за употребу
Bruksanvisning
Kullanım Talimatları
Інструкція для застосування
使用说明书
دستورالعمل‌های استفاده

REF: ONXI2
REF: ONXI3

CE
2797

 **On-X** life technologies, inc.™
designed for life

TABLE OF CONTENTS

INTRODUCTION	2
CUSTOMER SERVICE	2
INDICATIONS.....	2
CONTRAINDICATIONS	2
Sizers	2
Aortic Replica Sizers	2
Instrument Handle	3
Rotators.....	4
Mitral Rotator	4
Leaflet Probe	4
Instrument Tray.....	4
USE OF INSTRUMENTS	4
PREOPERATIVE AND OPERATIVE PRECAUTIONS AND WARNINGS	5
INSTRUMENT REPROCESSING	5
Reprocessing Limitation	5
Point of Use	5
Cleaning Preparation.....	5
Cleaning: Manual Method	6
Cleaning: Automated Method.....	6
Disinfection	6
Drying.....	6
Inspection	6
Packaging	7
Sterilization	7
Storage.....	7
RETURNED GOODS.....	8
DISCLAIMER OF WARRANTIES	8
SYMBOLS GLOSSARY DEFINITION	8

INTRODUCTION

This booklet is intended to aid the clinician by providing instructions for the handling and use of the ON-X LTI On-X® Prosthetic Heart Valve Instruments (Ref: ONXI2 and ONXI3). The instruments are comprised of sizers, rotators, instrument handle(s), and a leaflet motion probe(s). Upon receipt, if any instrument is damaged, DO NOT USE. Call ON-X LTI Customer Service to arrange for a replacement.

CUSTOMER SERVICE

ON-X LTI Customer Service representatives can be reached at (888) 339-8039 (USA and Canada only), (512) 339-8039 or email address ONX@ONXLTI.COM

INDICATIONS

The Instruments are intended to be used to facilitate the implantation of On-X® Prosthetic Heart Valves and On-X® Ascending Aortic Prostheses only.

CONTRAINDICATIONS

The Instruments should not be used for any purpose other than the implantation of On-X® Prosthetic Heart Valves or On-X® Ascending Aortic Prostheses.

DESCRIPTION AND FUNCTIONS OF THE INSTRUMENTS

The instrument kit contains a sizer and a rotator for each size of the On-X® Prosthetic Heart Valve and On-X® Ascending Aortic Prosthesis. Aortic replica sizers are included for size 19 through 27/29 mm aortic valves. A universal instrument handle and leaflet probe are also included. The instruments are reusable.

Sizers

The sizers are used to gauge the resulting tissue annulus diameter after the annulus is prepared for implant. The sizers are designed to be used for both aortic and mitral valve sizing and include size designations. The sizers are cylindrical for size 19 through 25 mm valves (Figure 1A) and conical for size 27/29 and 31/33 mm valves (Figure 1B). To facilitate size selection, refer to Table 1.

The sizer has a bendable stem on each end. The sizer has a bendable stem on each end. Repeated bending of the stems may cause metal fatigue. The effects of metal fatigue can be seen by the user in the form of readily apparent surface cracks on the stem. To avoid instrument fracture during use, the stem should be inspected for surface cracks before and after each time it is bent. If metal fatigue surface cracks are present, the sizer should be discarded and replaced. Contact ON-X LTI Customer Service to receive replacement sizers.

Aortic Replica Sizers

The aortic replica sizer (Figure 1A) models the On-X® aortic valve profile. It is used after sizing to assure fit of the aortic valve without obstruction of the coronary arteries. Aortic replica sizers are provided for size 19 through 25 mm aortic valves, where the valve sewing cuff is intended to remain supra-annular. The size 27/29 mm aortic valve is an intra-annular design and its replica sizer is designed to model this positioning.


Table 1

Size (mm)	Valve/Type	Sizer Choice	
		Sizer Type	Use Replica Sizer
19	Aortic*	Cylindrical	YES
21	Aortic*	Cylindrical	YES
23	Aortic*	Cylindrical	YES
25	Aortic*	Cylindrical	YES
27/29	Aortic*	Conical	YES
23	Mitral	Cylindrical	NO
25	Mitral	Cylindrical	NO
27/29	Mitral	Conical	NO
31/33	Mitral	Conical	NO
25/33	Mitral Conform - X	Cylindrical or Conical	NO

*For On-X® Prosthetic Heart Valves and On-X® Ascending Aortic Prostheses

ONX12 Sizers are made with transparent "smoke" colored plastic and are designed to be the same size as the corresponding heart valves. ONX13 Sizers are made with a transparent "gold" colored plastic and are slightly oversized to estimate annular size with pledgeted sutures in place.

Instrument Handle

The instrument handle (Figure 2) facilitates holding the valve or the rotator during surgery. The instrument handle is comprised of a grip, a bendable stem and a tip.

The instrument handle tip is inserted into the valve holder while the valve is still in the package inner container. The tip is inserted into the valve holder by placing it directly into the slot on top of the valve holder. It snaps into place after the application of a light insertion force. Upon snapping into position, the valve and the holder are firmly retained by the instrument handle. Removal of the valve from the inner container is performed after the valve holder is snapped onto the instrument handle.


Figure 2: Instrument Handle

The instrument handle has a bendable stem. Repeated bending of the stem may cause metal fatigue. The effects of metal fatigue can be seen by the user in the form of readily apparent surface cracks on the stem. To avoid instrument fracture during use, the stem should be inspected for surface cracks before and after each time it is bent. If metal fatigue surface cracks are present, the instrument handle should be discarded and replaced. Contact ON-X LTI Customer Service to receive a replacement instrument handle.

Rotators

The valve rotator (Figure 3) is used for reorienting an in situ valve and may be used to verify leaflet mobility. The valve rotator consists of a plastic head with a centrally located leaflet mobility probe and an attached handle. The valve rotator is properly oriented for insertion into the valve when the cross-bar on the head is aligned with the leaflet pivot axis and the probe tip is inserted into the central orifice between the leaflets.


Figure 3: Valve Rotator

Mitral Rotator

The mitral rotator (Figure 4) is used for reorienting an in situ valve in the mitral position only. The mitral rotator is distinguishable from the valve rotator by its shorter body and lack of a leaflet probe.

The rotators may be used with or without the instrument handle attached. To attach a rotator to the instrument handle, insert the instrument handle tip directly into the slot on the end of the rotator handle. The rotator snaps into place after the application of a light insertion force. The rotator must be removed from the handle after use prior to cleaning. A force greater than the insertion force is required to remove the rotator from the instrument handle.


Figure 4: Mitral Rotator

Leaflet Probe

The leaflet probe (Figure 5) is a flexible rod with tapered ends. The leaflet probe may be used to gently move the leaflets to verify that they open and close freely.


Figure 5: Leaflet Probe

Instrument Tray

The instrument tray is optional and may be used for storage and transportation of On-X® Valve instruments. The tray can withstand standard and flash steam sterilization. The instruments may be sterilized while in this tray. **WARNING: The tray, by itself, does not maintain sterility of the instruments after sterilization. To maintain sterility of instruments, the tray with instruments must be sterilized in conjunction with a wrapping material that allows steam sterilization and maintains sterility.**

USE OF INSTRUMENTS

Additional instructions for the use of each instrument are included in the Instructions for Use for the Prosthetic Heart Valve and Ascending Aortic Prosthesis. Those instructions must also be followed.

The instruments are reusable and must be separately cleaned and sterilized prior to each use.

The life cycle for the On-X instruments is 10 years from the date of manufacture or 100 reprocessing cycles, whichever comes first. Instrument usage and reprocessing should be monitored in accordance with the healthcare facility's record keeping protocol for equipment maintenance.

PREOPERATIVE AND OPERATIVE

PRECAUTIONS AND WARNINGS

CAUTION: Use these instruments to facilitate the implantation of the On-X® Prosthetic Heart Valves or On-X® Ascending Aortic Prostheses only.

CAUTION: Test the leaflet mobility only with the On-X® leaflet probe or the leaflet probe on the end of the On-X® rotator.

CAUTION: Use only the On-X® rotator to rotate the valve in situ. Do not attempt to rotate the valve if any significant resistance to rotation is encountered.

WARNING: Do not contact the carbon components of the valve with any metallic or abrasive instruments.

WARNING: Under no circumstances should leveraged force be applied to the valve housing or leaflets, as serious damage may occur.

INSTRUMENT REPROCESSING

This section specifies reprocessing instructions describing the methods to thoroughly clean and steam sterilize the instruments.

WARNING: Chemical disinfectants that contain phenolics, acetone, ammonium compounds, iodophors and formaldehyde should not be used to disinfect the instruments.

WARNING: Do not expose instruments with plastic components to temperatures at or above 350°F (177°C), otherwise damage to plastic component may occur.

WARNING: Do not use or repair a damaged instrument. Contact ON-X LTI customer service to receive a replacement instrument.

Reprocessing Limitation

The life cycle for the On-X instruments is 10 years from the date of manufacture or 100 reprocessing cycles, whichever comes first. Instrument usage and reprocessing should be monitored in accordance with the healthcare facility's record keeping protocol for equipment maintenance.

The plastic valve holder that comes pre-mounted to the valve is a single use device and should be discarded after valve implantation.

Metallic parts of the instruments are made of titanium and the plastic instruments are made of polyphenylsulfone (common brand name is RADEL® R a product of Solvay Advanced Polymers).

Point of Use

Within a maximum of two (2) hours after use, remove excess soil on instrument with a disposable cloth/paper wipe. Soiled instruments should not be allowed to dry, as dried contaminants are more difficult to remove during cleaning. If transfer is likely to be delayed, place instruments in covered container with appropriate detergent (e.g. Enzol® Enzymatic Detergent or equivalent) to delay drying.

Cleaning Preparation

It is recommended that the instruments be cleaned as soon as reasonably practical after use. Before cleaning, remove the valve holder or rotator from the instrument handle. Discard the valve holder.

Accessories that may be used during the cleaning preparation and cleaning process are surgical scrub type brushes, chenille pipe cleaners, cotton tipped applicators and cloth wipes.

Preparation of instruments for cleaning by soaking or ultrasonic cleaning is recommended. The use of a brush to remove gross contamination is also recommended.

If instruments are heavily soiled, immerse in a container of an appropriate cleaner (e.g. Enzol Enzymatic Detergent prepared with one (1) ounce to one (1) gallon tap water). Using a soft bristle brush, clean the instrument thoroughly, removing all visible soil. Enzyme products will not degrade the instruments. Rinse in warm running tap water until all traces of cleaning solution are removed. Repeat if necessary to remove all visible soil.

Cleaning: Manual Method

Wash all surfaces of the instrument with a prepared neutral, non enzymatic detergent solution (for example, Luminox at the manufacturers specified concentration) and a soft bristled brush until all visible traces of debris are removed.

Rinse instrument with water to remove all visible traces of detergent residue, rinse blind holes repeatedly.

Visually inspect instrument for any remaining soil or detergent residue. If found repeat cleaning procedure using fresh solution.

Cleaning: Automated Method

Place instruments in an automated washer (e.g. Steris 444 Washer/Disinfector or equivalent). Load instruments such that contact is avoided. For reference, a validated set of cleaning process parameters is given below:

Treatment	Time (mm:ss)	Temperature (°C)	Cleaning Agent
Enzymatic Wash	04:00	60	Enzol Enzyme Detergent (1 oz/gal)
Wash	02:00	60	Steris® Prolystica® 2X Concentrate Neutral (1/8 oz/gal)
Rinse	02:00	70	Not Applicable
Dry	15:00	80	Not Applicable

Visually inspect each instrument for remaining soil or moisture. If soil remains, repeat the cleaning process. If instruments are wet, use filtered compressed air or a lint-free wipe to dry.

WARNING: Do not expose instruments with plastic components to temperatures at or above 350°F (177°C), otherwise damage to plastic components may occur.

Disinfection

Use of chemical disinfectants is not required, although permissible. If used, chemical disinfectants such as alcohol, hydrogen peroxide, chlorine and glutaraldehyde will not damage or degrade the instruments. Disinfectants other than alcohol must be thoroughly rinsed from the surface of the instruments after use.

Disinfectants that contain phenolic, ammonium compounds, acetone, iodophors or formaldehyde should not be used to disinfect the instruments.

Drying

Do not expose instruments with plastic components to temperatures at or above 350°F (177°C), otherwise damage to plastic components may occur.

Inspection

Visually inspect instruments for damage and excessive wear. Sizers and instrument handles have bendable stems. Repeated bending of the stems may cause metal fatigue which can be seen in the form of readily apparent surface cracks on the stem. Instruments, either metal or plastic, exhibiting cracks, nicks or other damage should be replaced.

Do not use or repair a damaged or worn instrument. Contact ON-X LTI customer service to receive a replacement instrument.

Packaging

If an optional Instrument Tray is used, instruments should be placed into the identified locations and the lid attached. Double wrap instrument tray in CSR hospital wraps (such as Bioshield II 40" x 60") using an envelope or square wrap technique and secured with autoclave tape. If an Instrument Tray is not used, instruments should be stored in a manner that prevents bending or crushing.

CAUTION: The use of paper/plastic pouches is not recommended for the packaging of instruments, as condensate may remain trapped in the pouch.

Sterilization

ONXI2 and ONXI3 instruments are considered "Critical Items" that may be introduced directly into the blood stream or other normally sterile areas of the body. For sterilization, trays should be wrapped with two layers of CSR hospital wraps and secured with autoclave tape. They must be steam sterilized using the validated minimum time and temperature parameters specified below:

Validated Sterilization Cycle Parameters	
Method:	Steam
Cycle Type:	Pre-vacuum
Temperature:	132° C (270° F)
Exposure Time:	3 minutes
Dry Time:	20 minutes

These parameters are only valid with sterilization equipment that is properly maintained and calibrated. The water quality required for sterilization shall meet recommendations provided by the sterilization equipment manufacturer.

WARNING: Do not expose instruments with plastic components to temperatures at or above 350°F (177°C), otherwise damage to plastic component may occur.

Storage

If the instruments are not to be used immediately after sterilization, they may be sealed within a plastic dust cover (maintenance cover) and then stored in a clean, cool and dry environment. Heat sealing or a similarly effective method may be used to seal the maintenance cover.

The shelf life of the package is event-related and depends on the wrapping material and the manner in which it is stored. Therefore, the instruments should be handled in a manner that will prevent the sterile barrier from getting wet, being punctured or becoming compromised in some other manner. If the sterile barrier has been compromised, then the instruments must be re-cleaned and re-sterilized.

If a tray is not available, the instruments should be stored in a manner that will prevent them from becoming bent or crushed.

Warning: any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

RETURNED GOODS


Prior authorization from Customer Service is required for the return of any product. For any questions regarding the instruments or for the return authorization, please contact ON-X LTI Customer Service.

CAUTION: Federal (U.S.) law restricts this device to sale by or on the order of a physician.

DISCLAIMER OF WARRANTIES

ON-X LTI warrants only that the product shall conform to ON-X LTI's standard specifications. No other warranty is made by ON-X LTI concerning the function of the product in use, and ON-X LTI assumes no risk whatsoever as to the results of the use of the product. The entire risk with use of the product is that of the buyer. ON-X LTI disclaims all other warranties respecting the product, expressed or implied, including but not limited to those related to the product's merchantability or fitness for a particular purpose. ON-X LTI shall not be liable for any direct, special, consequential or incidental loss, damage or expense related to the use of the product. No person has any authority to alter any of these conditions or to bind ON-X LTI to any additional responsibility or warranty in connection with the use of the product.

SYMBOLS GLOSSARY DEFINITION

SYMBOL	SYMBOL TITLE	DEFINITION	STANDARD TITLE AND REFERENCE
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1 Medical Devices – Symbols To Be Used with Medical Device Labels, Labeling, and Information to be Supplied – Part 1: General requirements Reference no. 5.1.1
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1 Medical Devices – Symbols To Be Used with Medical Device Labels, Labeling, and Information to be Supplied – Part 1: General requirements Reference no. 5.1.5
	Consult instructions for use	Indicates the need for the user to consult the instructions for use	ISO 15223-1 Medical Devices – Symbols To Be Used with Medical Device Labels, Labeling, and Information to be Supplied – Part 1: General requirements Reference no. 5.4.3
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified	ISO 15223-1 Medical Devices – Symbols To Be Used with Medical Device Labels, Labeling, and Information to be Supplied – Part 1: General requirements Reference no. 5.1.6
	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.	ISO 15223-1 Medical Devices – Symbols To Be Used with Medical Device Labels, Labeling, and Information to be Supplied – Part 1: General requirements Reference no. 5.2.8
	Authorized representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1 Medical Devices – Symbols To Be Used with Medical Device Labels, Labeling, and Information to be Supplied – Part 1: General requirements Reference no. 5.1.2
	Date of Manufacture	Indicates the date that the medical device was manufactured	ISO 15223-1 Medical Devices – Symbols To Be Used with Medical Device Labels, Labeling, and Information to be Supplied – Part 1: General requirements Reference no. 5.1.3
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	ISO 15223-1 Medical Devices – Symbols To Be Used with Medical Device Labels, Labeling, and Information to be Supplied – Part 1: General requirements Reference no. 5.4.4

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On-X Life Technologies, Inc.
1300 East Anderson Lane, Bldg. B, Austin, Texas 78752 U.S.A.
JOTEC GmbH, a fully owned subsidiary of CryoLife, Inc.